

Routine Monitoring Seal-Check

4a medical[®]
supporting healthcare




Routine Monitoring Seal-Check

4a medical[®]
supporting healthcare



Acceptable
sealing for the
specified sealing width
in accordance with EN
ISO 11607-2.



6.061.017

PEEL ↑

| Sealing Machine · Siegelgerät Apareil de Scellage · Máquina selladora | Location · Standort Site · Ubicación | For Approval · Freigabe Autorisation · Habilitación | Date · Datum Date · Fecha | Remarks · Bemerkungen Remarques · Nota |
|--|---|--|------------------------------|---|
| | | | | |

6.061.017


hawa

Seal Check

T +49 (0) 6261/9770-0 · F +49 (0) 6261/62015 · info@hawa.com · www.hawa.com

No usare se il papel o el sellado se encuentra deteriorado
Non utilizzare se la confezione è danneggiata

Channels or open sealings.



PEEL ↑

| Sealing Machine · Siegelgerät Apareil de Scellage · Máquina selladora | Location · Standort Site · Ubicación | For Approval · Freigabe Autorisation · Habilitación | Date · Datum Date · Fecha | Remarks · Bemerkungen Remarques · Nota |
|--|---|--|------------------------------|---|
| 6.061.017 | | | | |


▲▲ Seal Check

hawo
T +49 (0) 6261/9770-0 · F +49 (0) 6261/62015 · info@hawo.com · www.hawo.com

▲▲

No usare se il papel o el sellado se encuentra deteriorado
Non utilizzare se la confezione è danneggiata

Delamination or separation of material/temperature or contact pressure too low.



6.061.017

PEEL ↑

hawo

| Sealing Machine · Siegelgerät Appareil de Scellage · Máquina selladora | Location · Standort Site · Ubicación | For Approval · Freigabe Autorisation · Habilitación | Date · Datum Date · Fecha | Remarks · Bemerkungen Remarques · Nota |
|---|---|--|------------------------------|---|
| | | | | |


6.061.017

Seal Check

hawo

T +49 (0) 6261/9770-0 · F +49 (0) 6261/62015 · info@hawo.com · www.hawo.com

No usare el papel o el sellado se encuentra deteriorado
Non utilizzare se la confezione è danneggiata



punctures
or tears

6.061.017

PEEL ↑

| | | | | |
|---|---|--|------------------------------|---|
| Sealing Machine · Siegelgerät Appareil de Scellage · Máquina selladora | Location · Standort Site · Ubicación | For Approval · Freigabe Autorisation · Habilitación | Date · Datum Date · Fecha | Remarks · Bemerkungen Remarques · Nota |
|---|---|--|------------------------------|---|


Seal Check

hawo

T +49 (0) 6261/9770-0 · F +49 (0) 6261/62015 · info@hawo.com · www.hawo.com

No usar si el papel o el sellado se encuentra deteriorado
Non utilizzare se la confezione è danneggiata

Temperature or contact pressure is too high.



PEEL ↑

| Sealing Machine · Siegelgerät Apareil de Scellage · Máquina selladora | Location · Standort Site · Ubicación | For Approval · Freigabe Autorisation · Habilitación | Date · Datum Date · Fecha | Remarks · Bemerkungen Remarques · Nota |
|--|---|--|------------------------------|---|
| 6.061.017 | | | | |